**Hager & Meisinger GmbH**

Hansemannstr. 10  
41468 Neuss • Germany  
Tel.: +49 (0) 21 31-20 120  
Fax: +49 (0) 21 31-20 12 222  
Internet: [www.meisinger.de](http://www.meisinger.de)  
E-mail: [info@meisinger.de](mailto:info@meisinger.de)

**Meisinger USA, L.L.C.**

12276 San Jose Boulevard, Suite 110  
Jacksonville, Florida 32223 • USA  
Tel.: +1 (904) 260-5040  
Fax: +1 (904) 260-5212  
Internet: [www.meisingerusa.com](http://www.meisingerusa.com)  
E-mail: [info@meisingerusa.com](mailto:info@meisingerusa.com)

HMV97-03-1/2/3-25

**Split-Control**  
**Bone Spreading**  
**and**  
**Condensing System**

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In the case of a horizontal alveolar ridge resorption, Bone-Spreading offers the possibility of gentle dilatation of the existing bone substance and may be deployed prior to all current dental implant systems.

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The preliminary deployment of Split-Control also perfectly prepares the bone cavity in the case of a horizontally resorbed alveolar ridge via dilatation.

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
Cortical osteotomy with an separating disc to promote bone-spreading.

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The pilot bur produces a subdimensional bone cavity (take notice of diameter and depth of the drilling).

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The deployment of a series of threadformers with increasing diameters leads to bonespreading and condensation (order of diameter and codification of the instruments must be respected). With the help of a carrier (short/long) the threadformers may be screwed in cautiously, and if necessary, with the help of the "driver" (pay attention to depth of insertion).

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The precondensed implant cavity increases the primary stability irrespective of the subsequently deployed implant system.

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The particular properties of Split-Control enable a substance-saving densification of cancellous bone of sufficient horizontal dimension and therefore prepare the implant cavity perfectly for the subsequent process of implant insertion.

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Split-Control can prepare the implant cavity perfectly by spreading the alveolar ridge and condensing the cancellous bone.

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In a narrow alveolar ridge situation this initial bur prevents the pilot bur from slipping.

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The pilot bur creates a sub-dimensional pilot drilling (take notice of diameter and depth of drilling).

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A sub-dimensional bone cavity for bone densities D3/4 is created by the conical expansion burs. For bone condensing purposes of implant cavities of more than 2.7 mm diameter, use the expansion bur A1 and the threadformer A1. For implant cavities exceeding 3.3 mm diameter, the expansion bur D1 and the threadformer D1 must be

used. After deployment of the primary threadformer (A1 or D1) the following threadformers in order can be used without prior drilling.

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In non-atrophied bone, the use of the threadformers leads to condensation of the bone. Horizontally resorbed alveolar ridge is spread additionally. With the help of a carrier (short/long), the threadformers may be screwed in cautiously and if necessary, with the help of the "driver" (pay attention to depth of insertion).

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The precondensed (and if necessary dilated) implant cavity thus increases the primary stability irrespective of the type of implant to be inserted subsequently.

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Art.-Nr.: CBS01

- Gentle bone-spreading through controlled deployment of non-ablative thread-formers
- Excellent alternative to the bone traumatism created by hammer and chisel preparation of the bone. Split-Control, a set of high definition screw-like tools, achieves a controlled and standardized bone-dilatation and condensation.
- Suitable for all bone resorption classes (1–5)
- High primary stability through bone-condensing
- Suitable for preliminary implant cavity preparation of all current dental implant systems

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